

# HEAD LUX

## FIBEROPTIC & LED OPTIONS

The extended family of Ecleris Headlights,  
Everlux & Headlux Products

Premium quality dual adjustment head band  
Optimized illumination



FIBEROPTIC HEADLUX

LED HEADLUX



# HEAD LUX

## FIBEROPTIC & LED OPTIONS

The extended family of Ecleris Headlights, Everlux & Headlux Products

High light transmission - Adjustable focus with extremely precise beam  
Light weight – adjustable – comfortable

### • FIBEROPTIC HEADLUX



Standard fiber optic cable connection (Storz type)

- Cold light
- Best combination with Ecleris EVERLUX Mini



### • LED HEADLUX



Battery operated:  
7 hours continuous use at maximum intensity



Adjustable focus	Min.	Max.
Diameter of light spot at 340 mm distance	19mm	70mm
Diameter of light spot at 420 mm distance	21mm	90mm
Illumination at 200 mm distance	≥ 35,000 lx	
Electrical specifications		
Illumination source	LED 20,000 hours self-life	
Rechargeable Battery	7.4V Li-ion polymer battery	
Battery Run-time	≥ 7 hours	
Charging time	≤ 3.5 hours	